



# KAPA Library Preparation Kits

## with Real-Time PCR Library Amplification

Illumina series

### Product Description

The KAPA Library Preparation Kit provides all of the enzymes and reaction buffers required for constructing libraries from fragmented dsDNA via the following steps:

- 1. End repair:** Produce blunt-ended, 5'-phosphorylated fragments.
- 2. A-tailing:** Add dAMP to the 3'-ends of the dsDNA library fragments.
- 3. Adaptor ligation:** Ligate dsDNA adaptors with 3'-dTMP overhangs to library fragments.
- 4. Library amplification:** Real-time amplification of library fragments carrying appropriate adaptor sequences on both ends.

Reaction buffers are supplied in convenient, concentrated "master mix" formats comprising all of the required reaction components except oligonucleotide adaptors and PCR primers. Similarly, a single enzyme mixture is provided for each step of the library construction process, reducing the number of pipetting steps.

High fidelity PCR is used to selectively enrich library fragments carrying appropriate adaptor sequences and to amplify the amount of DNA prior to sequencing. During PCR enrichment of libraries, minimizing amplification bias is critical to ensure more uniform sequence coverage. Amplification bias occurs when a DNA polymerase is unable to amplify all targets within a complex population of library DNA with equal efficiency. Bias is further exacerbated when libraries are over-amplified.

KAPA HiFi Real-Time PCR Library Amplification Kits are designed to address both sources of PCR-induced bias. The novel KAPA HiFi DNA Polymerase, engineered for high fidelity and processivity, is capable of balanced amplification of complex library DNA. Real-time monitoring of library amplification provides additional information required to optimize the number of amplification cycles and minimize over-amplification. The benefits of performing high fidelity, real-time PCR for next-generation sequencing library amplification include:

- Real-time monitoring of amplification allows precise control over the optimal number of PCR cycles.
- Real-time amplification workflows are amenable to automation.
- Real-time amplification plots provide quality metrics for individual enriched libraries, eliminating expensive and time-consuming post-enrichment gel electrophoresis and identifying inconsistencies in library preparation.
- Seamless integration with KAPA Library Quantification Kits.

KAPA HiFi Real-Time PCR Library Amplification Kits contain KAPA HiFi HotStart Real-Time PCR Master Mix (2X), a ready-to-use cocktail containing all components for PCR, except primers and template. The 2X Master Mix contains KAPA HiFi HotStart DNA Polymerase in a proprietary reaction buffer, dNTPs, MgCl<sub>2</sub> (2.5 mM at 1X), SYBR® Green I dye and stabilizers. Four fluorescent standards are supplied, and are used to define a window for optimal amplification (Figures 1 and 2).

This kit is primarily intended for the construction of genomic, paired-end, and paired-end multiplex (indexed/barcoded) Illumina DNA libraries, but may be used for other applications requiring efficient end-repair, A-tailing, ligation, and/or library amplification steps.

### Kit codes and components

<b>KK8220</b> 10 reactions	KAPA Library Preparation Kit with Real-Time PCR Library Amplification - Illumina <ul style="list-style-type: none"> <li>- End Repair Enzyme Mix (50 µL)</li> <li>- 10X End Repair Buffer with dNTPs (100 µL)</li> <li>- A-Tailing Enzyme (30 µL)</li> <li>- 10X A-Tailing Buffer (50 µL)</li> <li>- DNA Ligase (50 µL)</li> <li>- 5X Ligation Buffer (100 µL)</li> <li>- 2X KAPA HiFi HotStart Real-Time PCR Master Mix (250 µL)</li> <li>- 4 x Fluorescent Standards (1500 µL each)</li> </ul>
<b>KK8221</b> 50 reactions	KAPA Library Preparation Kit with Real-Time PCR Library Amplification - Illumina <ul style="list-style-type: none"> <li>- End Repair Enzyme Mix (250 µL)</li> <li>- 10X End Repair Buffer with dNTPs (500 µL)</li> <li>- A-Tailing Enzyme (150 µL)</li> <li>- 10X A-Tailing Buffer (250 µL)</li> <li>- DNA Ligase (250 µL)</li> <li>- 5X Ligation Buffer (500 µL)</li> <li>- 2X KAPA HiFi HotStart Real-Time PCR Master Mix (1250 µL)</li> <li>- 4 x Fluorescent Standards (1500 µL each)</li> </ul>

### Storage and handling

Store all components at -20 °C. KAPA HiFi HotStart Real-Time PCR Master Mix (2X) and fluorescent standards 1 – 4 are light sensitive and should be protected from light during storage, thawing, and reaction setup. Please refer to Section 6 for full details.

### Quick Notes

- Reaction components should be **mixed fresh** and used on the same day.
- Reaction cleanup between successive enzymatic reactions may be accomplished via your method of choice. We recommend Agencourt Ampure XP Beads or Qiagen MinElute Reaction Cleanup Kits.
- Optimal real-time library amplification corresponds to the region between fluorescent standard 1 and 3. The termination cycle number should be adjusted accordingly without the requirement for performing gel electrophoresis (see Figure 1, page 9).
- It is critical that the correct data acquisition temperature is adhered to, to minimize background fluorescence due to inter- and intra-primer interaction (see Table 1, page 12).
- When using primers that differ in sequence from those listed in Table 1, we recommend performing gradient PCR to optimize the annealing temperature.



## 1. End Repair

### Materials required but not supplied in this kit:

- Reaction tubes. Reactions may be assembled and processed in PCR plates, PCR tubes, or microcentrifuge tubes.
- Pipette tips. We strongly recommend the use of high quality filter-plugged tips to prevent contamination of reagents and library samples.
- Reaction cleanup. We recommend Agencourt AMPure XP Beads (cat. # A63881) or Qiagen MinElute Reaction Cleanup Kits (cat. # 28204) for this purpose. AMPure XP Beads require the preparation of an elution buffer (10 mM TRIS-Acetate, pH 8.0, reagent grade water).

#### 1.1 Assemble the End Repair reaction:

Water to 100 $\mu$ L	x $\mu$ L
10X End Repair Buffer	10 $\mu$ L
End Repair Enzyme Mix	5 $\mu$ L
1 - 5 $\mu$ g sheared dsDNA	1 - 85 $\mu$ L
<hr/>	
Total	100 $\mu$ L

#### 1.2 Incubate for 30 min @ 20 °C

#### 1.3 Proceed immediately to cleanup.

## End Repair Cleanup

We recommend either Agencourt AMPure XP Beads, **or** Qiagen MinElute Reaction Cleanup Kits. Below are suggested protocols for the recommended methods:

### AMPure XP Beads

#### 1.4 Ensure that the AMPure XP Beads are equilibrated to room temperature, and that they are thoroughly resuspended.

#### 1.5 Add AMPure Beads to the End Repair reaction:

End Repair reaction	100 $\mu$ L
<b>AMPure XP Beads</b>	<b>160 <math>\mu</math>L</b>
<hr/>	
Total	260 $\mu$ L

#### 1.6 Mix thoroughly on a vortex mixer or by pipetting up and down at least ten times.

#### 1.7 Incubate at room temperature for 15 minutes to allow DNA to bind to the beads.

#### 1.8 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

#### 1.9 Carefully remove and discard 255 $\mu$ L of the liquid. Take care not to disturb or discard any of the beads. Some liquid may remain visible in the tube/well.

#### 1.10 Keeping the tube/plate on the magnetic stand and without disturbing the beads, wash the beads in 200 $\mu$ L of 80% EtOH for at least 30 seconds.

#### 1.11 Carefully remove and discard the ethanol without disturbing the beads, and repeat the process for a total of 2 washes in 80% EtOH.

#### 1.12 Remove the tube/plate from the magnetic stand, and allow the beads to dry at room temperature for 15 minutes.

#### 1.13 Resuspend the beads thoroughly in 32.5 $\mu$ L elution buffer, and incubate at room temperature for 2 minutes to release the DNA from the beads.

#### 1.14 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

#### 1.15 Recover the DNA in 30 $\mu$ L of supernatant and transfer to the tube/well in which you intend to perform the A-tailing reaction.



## 1. End Repair (cont.)

OR

### Qiagen MinElute Reaction Cleanup Kit

The procedure for reaction cleanup using Qiagen MinElute columns described below is abbreviated and is intended for users who are already familiar with the method. Please consult the documentation provided by the supplier for full details.

- 1.4 If the End Repair reaction was performed in a PCR tube/plate, transfer the End Repair reaction to a microcentrifuge tube before proceeding.
- 1.5 Follow the MinElute protocol:

End Repair reaction	100 $\mu$ L
<b>Add Buffer ERC</b>	<b>300 <math>\mu</math>L</b>
<hr/>	
Total	400 $\mu$ L

- 1.6 Apply the mixture to a column and centrifuge or apply vacuum as appropriate. Discard the flow-through.
- 1.7 Wash with 750  $\mu$ L Buffer PE. Discard the flow-through.
- 1.8 Centrifuge for 2 minute at  $\geq 10,000 \times g$  to remove all traces of ethanol.
- 1.9 To elute, transfer the column to a clean, sterile microcentrifuge tube and add 31  $\mu$ L buffer EB. Incubate for 1 minute at room temperature, and centrifuge to recover  $\sim 30 \mu$ L.

### **\*\*Safe Stopping Point\*\***

If you are not proceeding to A-Tailing immediately, the protocol can be safely stopped here. Store at  $-20 \text{ }^{\circ}\text{C}$  for up to seven days.



## 2. A-Tailing

2.1 Assemble the A-Tailing reaction:

Water	12 $\mu$ L
10X A-Tailing Buffer	5 $\mu$ L
A-Tailing Enzyme	3 $\mu$ L
End repaired DNA	30 $\mu$ L
<hr/>	
Total	50 $\mu$ L

2.2 Incubate for 30 min @ 30 °C

2.3 Proceed immediately to cleanup.

### A-Tailing Cleanup

We recommend either Agencourt AMPure XP Beads, **or** Qiagen MinElute Reaction Cleanup Kits. Below are suggested protocols for the recommended methods:

#### AMPure XP Beads

2.4 Ensure that the AMPure XP Beads are equilibrated to room temperature, and that they are thoroughly resuspended.

2.5 Add AMPure XP Beads to the End Repair reaction:

A-Tailing reaction	50 $\mu$ L
<b>Add AMPure XP Beads</b>	<b>90 <math>\mu</math>L</b>
<hr/>	
Total	140 $\mu$ L

2.6 Mix thoroughly on a vortex mixer or by pipetting up and down at least ten times.

2.7 Incubate at room temperature for 15 minutes to allow DNA to bind to the beads.

2.8 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

2.9 Carefully remove and discard 135  $\mu$ L of the liquid. Take care not to disturb or discard any of the beads. Some liquid may remain visible in the tube/well.

2.10 Keeping the tube/plate on the magnetic stand and without disturbing the beads, wash the beads in 200  $\mu$  L of 80% EtOH for at least 30 seconds.

2.11 Carefully remove and discard the ethanol without disturbing the beads, and repeat the process for a total of 2 washes in 80% EtOH.

2.12 Remove the tube/plate from the magnetic stand, and allow the beads to dry at room temperature for 15 minutes.

2.13 Resuspend the beads thoroughly in 32.5  $\mu$ L elution buffer, and incubate at room temperature for 2 minutes to release the DNA from the beads.

2.14 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

2.15 Recover the DNA in 30  $\mu$ L of supernatant and transfer to the tube/well in which you intend to perform the adaptor ligation reaction.



## 2. A-Tailing (cont.)

OR

### Qiagen MinElute Reaction Cleanup Kit

The procedure for reaction cleanup using Qiagen MinElute columns described below is abbreviated and is intended for users who are already familiar with the method. Please consult the documentation provided by the supplier for full details.

- 2.4 If the end repair reaction was performed in a PCR tube/plate, transfer the End Repair reaction to a microcentrifuge tube before proceeding.
- 2.5 Follow the MinElute protocol:

A-Tailing Reaction	50 $\mu$ L
<b>Add Buffer ERC</b>	<b>300 <math>\mu</math>L</b>
<hr/>	
Total	350 $\mu$ L

- 2.6 Apply the mixture to a column and centrifuge or apply vacuum as appropriate. Discard the flow-through.
- 2.7 Wash with 750  $\mu$ L Buffer PE. Discard the flow-through.
- 2.8 Centrifuge for 2 minute at  $\geq 10,000 \times g$  to remove all traces of ethanol.
- 2.9 To elute, transfer the column to a clean, sterile microcentrifuge tube and add 31  $\mu$ L buffer EB. Incubate for 1 minute at room temperature, and centrifuge to recover  $\sim 30 \mu$ L.

### **\*\*Safe Stopping Point\*\***

If you are not proceeding to Adaptor Ligation immediately, the protocol can be safely stopped here. Store at  $-20 \text{ }^{\circ}\text{C}$  for up to seven days.



### 3. Adaptor Ligation

3.1 Assemble the Adaptor Ligation reaction:

5X Ligation Buffer	10 $\mu$ L
DNA Ligase	5 $\mu$ L
DNA Adaptor (30 $\mu$ M)*	5 $\mu$ L
A-Tailed DNA	30 $\mu$ L
<hr/>	
Total	50 $\mu$ L

\*Not supplied

3.2 Incubate for 15 min @ 20 °C

3.3 Proceed immediately to cleanup.

#### Adaptor Ligation Cleanup

##### Note:

- If using AMPure XP Beads for DNA recovery and purification, we recommend that you follow the protocol outlined below to perform a total of **TWO** successive cleanup procedures after adaptor ligation. If you are using the MinElute Reaction Cleanup Kit, then a **single** cleanup procedure is recommended.

#### AMPure XP Beads

##### First AMPure XP Bead Cleanup

3.4 Ensure that the AMPure XP Beads are equilibrated to room temperature, and that they are thoroughly resuspended.

3.5 Add AMPure Beads to the End Repair reaction:

Adaptor Ligation reaction	50 $\mu$ L
<b>Add AMPure XP Beads</b>	<b>50 <math>\mu</math>L</b>
<hr/>	
Total	100 $\mu$ L

3.6 Mix thoroughly on a vortex mixer or by pipetting up and down at least ten times.

3.7 Incubate at room temperature for 15 minutes to allow DNA to bind to the beads.

3.8 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

3.9 Carefully remove and discard 95  $\mu$ L of the liquid. Take care not to disturb or discard any of the beads. Some liquid may remain visible in the tube/well.

3.10 Keeping the tube/plate on the magnetic stand and without disturbing the beads, wash the beads in 200  $\mu$ L of 80% EtOH for at least 30 seconds.

3.11 Carefully remove and discard the ethanol without disturbing the beads, and repeat the process for a total of 2 washes in 80% EtOH.

3.12 Remove the tube/plate from the magnetic stand, and allow the beads to dry at room temperature for 15 minutes.

3.13 Resuspend the beads thoroughly in 52.5  $\mu$ L elution buffer, and incubate at room temperature for 2 minutes to release the DNA from the beads.

3.14 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.

3.15 Recover the DNA in 50  $\mu$ L of supernatant and transfer to the tube/well in which you intend to perform the second cleanup procedure (see below).



### 3. Adaptor Ligation (cont.)

#### Second AMPure XP Bead Cleanup

- 3.16 Ensure that the AMPure XP Beads are equilibrated to room temperature, and that they are thoroughly resuspended.  
3.17 Add AMPure Beads to the End Repair reaction:

Library DNA from first cleanup	50 $\mu$ L
<b>Add AMPure XP Beads</b>	<b>50 <math>\mu</math>L</b>
<hr/>	
Total	100 $\mu$ L

- 3.18 Mix thoroughly on a vortex mixer or by pipetting up and down at least ten times.  
3.19 Incubate at room temperature for 15 minutes to allow DNA to bind to the beads.  
3.20 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.  
3.21 Carefully remove and discard 95  $\mu$ L of the liquid. Take care not to disturb or discard any of the beads. Some liquid may remain visible in the tube/well.  
3.22 Keeping the tube/plate on the magnetic stand and without disturbing the beads, wash the beads in 200  $\mu$ L of 80% EtOH for at least 30 seconds.  
3.23 Carefully remove and discard the ethanol without disturbing the beads, and repeat the process for a total of 2 washes in 80% EtOH.  
3.24 Remove the tube/plate from the magnetic stand, and allow the beads to dry at room temperature for 15 minutes.  
3.25 Resuspend the beads thoroughly in 32.5  $\mu$ L elution buffer, and incubate at room temperature for 2 minutes to release the DNA from the beads.  
3.26 Capture the beads by placing the tube/PCR plate on an appropriate magnetic stand at room temperature for 15 minutes or until the liquid is completely clear.  
3.27 Recover the DNA in 30  $\mu$ L of supernatant and transfer to the tube/well in which you intend to store the library DNA until you carry out size selection.

**OR**

#### Qiagen MinElute Reaction Cleanup Kit

The procedure for reaction cleanup using Qiagen MinElute columns described below is abbreviated and is intended for users who are already familiar with the method. Please consult the documentation provided by the supplier for full details.

- 3.4 If the end repair reaction was performed in a PCR tube/plate, transfer the End Repair reaction to a microcentrifuge tube before proceeding.  
3.5 Follow the MinElute protocol:

Adaptor Ligation reaction	50 $\mu$ L
<b>Add Buffer ERC</b>	<b>300 <math>\mu</math>L</b>
<hr/>	
Total	350 $\mu$ L

- 3.6 Apply the mixture to a column and centrifuge or apply vacuum as appropriate. Discard the flow-through.  
3.7 Wash with 750  $\mu$ L Buffer PE. Discard the flow-through.  
3.8 Centrifuge for 2 minute at  $\geq 10,000 \times g$  to remove all traces of ethanol.  
3.9 To elute, transfer the column to a clean, sterile microcentrifuge tube and add 31  $\mu$ L buffer EB. Incubate for 1 minute at room temperature, and centrifuge to recover  $\sim 30 \mu$ L.

#### **\*\*Safe Stopping Point\*\***

If you are not proceeding to Size Selection immediately, the protocol can be safely stopped here. Store at -20  $^{\circ}$ C for up to seven days.



## 4. Size Selection

It is important to remove unligated adaptor molecules prior to library amplification to prevent the formation of “adaptor dimers” and other short adaptor-derived molecules, which may cause problems downstream during cluster amplification and sequencing. For many sequencing libraries/protocols, it is also helpful to select a relatively narrow and precisely defined size range of library fragments. Depending on your needs and options, you may choose to perform this size selection via a variety of common methods including:


- “SPRI-bead” size selection.
- “Double-SPRI-bead” size selection.
- Manual agarose gel electrophoresis, excision, and purification.
- Automated DNA size selection and collection (e.g. Sage Science Pippin Prep™).

Due to the wide variety of viable alternatives for this procedure we do not provide any specific protocols for size selection.




## 5. Real-Time Library Amplification


### 5.1 Overview of real-time high fidelity amplification of next-generation DNA sequencing libraries.



**Figure 1.** Libraries are amplified using a SYBR® Green-based real-time, high fidelity PCR master mix (left panel). Four triplicate wells of the PCR plate contain fluorescent reference standards representing a range of distinct DNA concentrations. Reactions terminated between standards 1 and 3 represent the optimal library amplification range (grey box), depicted here from cycles 10 - 14. A gel image of a typical library stopped at different amplification cycles is shown in the right panel. Low- and high-molecular weight artifacts increase progressively with additional cycles.



**Figure 2a.** Superimposed amplification plots for reactions terminated at the lower bound (hashed line, cycle 10) or upper bound (solid line, cycle 14) of the targeted concentration range (grey box). Library amplification reactions should ideally be terminated anywhere within the indicated target concentration range.



**Figure 2b.** Example of real-time high fidelity amplification of multiple libraries. 20 libraries, spanning a ~64-fold concentration range (6 cycles), were simultaneously amplified and terminated after 14 cycles. 14 of the 20 libraries fall within the targeted amplification range. The remaining 6 libraries could either be used as is, noting that they may be outside the optimal concentration range, or they could be re-amplified individually or in high- or low-concentration groups.



## 5. Real-Time Library Amplification (cont.)

### 5.2 Preparation

**Note:** KAPA HiFi HotStart Real-Time PCR Master Mix (2X) and fluorescent standards 1 - 4 are light sensitive and should be protected from light.

- Thaw the primers required for PCR enrichment (PCR primers are not supplied in the KAPA Library Preparation Kit), a tube of KAPA HiFi HotStart Real-Time PCR Master Mix (2X), and fluorescent standards 1 - 4 at room temperature.
- Mix and briefly centrifuge the thawed KAPA HiFi HotStart Real-Time PCR Master Mix (2X), primer and fluorescent standards 1 - 4.
- Thaw, gently mix, and briefly centrifuge the adaptor-ligated, size-separated purified library DNA for 5 seconds at 600 xg.
- Pre-program the real-time thermal cycler with the appropriate cycling protocol (see Table 1) for each specific set of library amplification PCR primers.

### 5.3 Reaction setup

Each plate must contain a set of fluorescent standards 1 - 4 (50  $\mu$ L each loaded in triplicate) in addition to a single 50  $\mu$ L real-time PCR reaction for each library requiring amplification.

In order to maintain optimal library diversity it is necessary to add sufficient adaptor-ligated library DNA to each enrichment PCR reaction. The optimal cycle number is dependant on the volume and concentration of library material added to each 50  $\mu$ L PCR reaction. High background fluorescence may result if >100 ng dsDNA template is added per 50  $\mu$ L real-time PCR reaction.

Refer to Table 1 (page 12) for the suggested reaction setup for specific library preparation protocols.

Fluorescent standard set up requires the addition of 50  $\mu$ L of each fluorescent standard (1 - 4) in triplicate to wells of the real-time PCR plate.

Change tips after each pipetting step. Seal each reaction, mix gently and centrifuge briefly.

### 5.5 Cycling protocol

Refer to Table 1 for the thermal cycling protocol for specific library types.

- If conventional end-point PCR has previously been used successfully and the same amount and type of library is added to the KAPA HiFi HotStart Real-Time PCR reactions, then program the real-time thermocycler with the same number of cycles as previously used. If the optimal cycle number is not known refer to section 5.7 for more information.

### 5.6 Clean up PCR

- After enrichment PCR, clean up each reaction using either Agencourt AMPure XP beads (Beckman Coulter Genomics part # A63881) or Qiagen MinElute PCR Purification Kit (Qiagen, part # 28004).



## 5. Real-Time Library Amplification (cont.)

### 5.7 Validate library

Initially, the raw data (i.e., not background subtracted) linear real-time amplification plots can be used as a built-in quality metric to validate the level of amplification of each amplified library.

- If the linear amplification profile of the library is significantly below fluorescent standard 1 at the end of qPCR cycling, then it is unlikely that there will be sufficient library material to sequence after PCR purification.
- If the linear amplification profile of the library is significantly above fluorescent standard 3 at the end of qPCR cycling, then the library has been over-amplified. This may lead to 1) amplification bias, 2) higher error rates, and/or 3) the presence of chimeric PCR products.

This data is also useful as a quality control metric for identifying inconsistencies during library preparation between multiple libraries.

**NOTE:** The amplification plots can also be used in real-time to select the optimal cycle without a pre-programmed termination cycle. To do this:

1. Program 30 cycles into the real-time thermocycler.
2. After starting the real-time thermocycler, wait until the desired fluorescence of the library is achieved before terminating the real-time reaction.

**NOTE:** It is critical to terminate the reaction directly after data acquisition at 72°C **and before the tube ramps to 95°C** for the start of the next cycle. This will ensure that the enriched library DNA remains double-stranded for efficient downstream purification.

To verify the size of the PCR enriched fragments, check the size distribution by performing gel electrophoresis.

### 5.8 Library Quantification

Accurate quantification of amplifiable library molecules is critical for the efficient use of next-generation sequencing platforms. Overestimation of library concentration results in lower cluster density after bridge PCR. Underestimation of library concentration results in too many clusters on the flow cell, which can lead to poor cluster resolution. Both scenarios result in suboptimal sequencing capacity. Accurate library quantification is equally important when pooling indexed libraries for multiplexed sequencing to ensure equal representation of each library.

Integrate KAPA HiFi Real-Time PCR Library Amplification Kit with the appropriate KAPA Library Quantification Kit (KK4824, KK4835, KK4844, KK4854) to accurately quantify the number of PCR-competent molecules. If libraries have been terminated between fluorescent standards 1 – 3, a single 1:1,000 dilution of each library will be required for library quantification using the KAPA Library Quantification Kits.



## 5. Real-Time Library Amplification (cont.)

**Table 1. Recommended reaction setup and cycling parameters for KAPA HiFi HotStart Real-Time PCR Master Mix (2X) reactions:**

Library	Component	Final Conc.	Volume/50 µL rxn	Cycling Protocol
Genomic DNA	PCR grade water		As needed	Denaturation 45 sec at 98 °C
	2X KAPA HiFi HS RT-PCR MM	1X	25 µL	Cycling* 15 sec at 98 °C
ChIP	PCR Primer 1.1	500 nM	1 µL	30 sec at 65 °C
	PCR Primer 2.1	500 nM	1 µL	30 sec at 72 °C
	Library DNA		As needed	Final Extension 1 min at 72 °C
Paired-End	PCR grade water		As needed	Denaturation 45 sec at 98 °C
	2X KAPA HiFi HS RT-PCR MM	1X	25 µL	Cycling* 15 sec at 98 °C
	PE PCR Primer 1.0	500 nM	1 µL	30 sec at 65 °C
	PE PCR Primer 2.0	500 nM	1 µL	30 sec at 72 °C
	Library DNA		As needed	Final Extension 1 min at 72 °C
Paired-End Multiplex	PCR grade water		As needed	Denaturation 45 sec at 98 °C
	2X KAPA HiFi HS RT-PCR MM	1X	25 µL	Cycling* 15 sec at 98 °C
	PE PCR Primer InPE 1.0	500 nM	1 µL	30 sec at 65 °C
	PE PCR Primer InPE 2.0	10 nM	1 µL	30 sec at 72 °C
	PCR Primer Index 1 - 12	500 nM	1 µL	30 sec at 72 °C
	Library DNA		As needed	Final Extension 1 min at 72 °C
TruSeq™ DNA	2X KAPA HiFi HS RT-PCR MM	1X	25 µL	Denaturation 45 sec at 98 °C
	PCR Primer Cocktail (PPC)	500 nM each	5 µL	Cycling* 15 sec at 98 °C
	Library DNA		20 µL	30 sec at 60 °C 30 sec at 72 °C
				Final Extension 1 min at 72 °C

\* The optimal cycling number will depend upon the volume and concentration of adaptor-ligated, size separated, purified library DNA added to each enrichment PCR reaction. Typically this is in the 10-18 cycle range but may require optimization.



## Appendix A: Composition of Materials

### 10X KAPA End Repair Buffer - contains dNTPs and ATP

500 mM Tris-HCl  
100 mM MgCl<sub>2</sub>  
100 mM DTT  
10 mM ATP  
4 mM dATP  
4 mM dCTP  
4 mM dGTP  
4 mM dTTP  
pH 7.5 @ 25 °C

### KAPA End Repair Enzyme Mix

3,000 U/mL T4 DNA Polymerase  
10,000 U/mL T4 Polynucleotide Kinase  
Supplied in: 100 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM DTT, 0.1% Triton X-100, 50% Glycerol, pH7.4 @ 25 °C

### 10x KAPA A-Tailing Buffer - contains dATP

100 mM Tris-HCl  
100 mM MgCl<sub>2</sub>  
500 mM NaCl  
10 mM DTT  
2 mM dATP  
pH 7.9 @ 25 °C

### KAPA A-Tailing Enzyme

5,000 U/mL (0.5 mg/mL) Klenow Fragment (3' → 5' exo-)  
Supplied in: 20 mM Tris-HCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.5 @ 25 °C

### 5X KAPA Ligation Buffer

330 mM Tris-HCl  
50 mM MgCl<sub>2</sub>  
5 mM DTT  
5 mM ATP  
30% PEG 6000  
pH 7.6 @ 25 °C

### KAPA DNA Ligase

600,000 U/mL (2 mg/mL) T4 DNA Ligase  
Supplied in: 10 mM Tris-HCl, 50 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.5 @ 25 °C

### KAPA HiFi HotStart Real-Time PCR Master Mix

KAPA HiFi HotStart DNA Polymerase is an antibody-based hot start formulation of KAPA HiFi DNA Polymerase, a novel B-family DNA polymerase engineered for increased processivity and high fidelity. KAPA HiFi HotStart DNA Polymerase has 5' → 3' polymerase and 3' → 5' exonuclease (proofreading) activities, but no 5' → 3' exonuclease activity. The strong 3' → 5' exonuclease activity results in superior accuracy during DNA amplification. The error rate of KAPA HiFi HotStart DNA Polymerase is calculated at 1 error in 3.54 x 10<sup>6</sup> bases covered (2.82 x 10<sup>-7</sup>). The SYBR® Green I intercalating dye contained within the master mix does not affect the fidelity of the polymerase. DNA fragments generated with KAPA HiFi HotStart ReadyMix may be used for routine downstream analyses or applications, including restriction enzyme digestion and sequencing. PCR products generated with KAPA HiFi HotStart ReadyMix are blunt-ended, but may be 3'-dA-tailed for cloning into TA cloning vectors.



## 6. Storage, handling and specifications

### 6.1 Shipping, storage and handling

KAPA Library Preparation Kits are shipped on dry ice or ice packs, depending on the country of destination. Upon receipt, store the entire kit at -20 °C in a constant-temperature freezer. When stored under these conditions and handled correctly, all kit components will retain full activity until the expiry date indicated on the kit.

Please note that certain components in KAPA Library Preparation Kits (e.g. End Repair Mix, DNA Ligase, buffers containing dNTPs and/or ATP, etc.) are particularly sensitive to temperature and freeze-thaw cycles, and should be handled with special care. Minimize exposure of the KAPA HiFi HotStart Real-Time PCR Master Mix (2X) and Fluorescent Standards 1 - 4 to direct light. Exposure to direct light for an extended period of time may result in loss of fluorescent signal intensity.

The KAPA HiFi HotStart Real-Time PCR Master Mix (2X) contains isostabilizers and may not freeze solidly, even when stored at -20 °C. Nevertheless, always ensure that the KAPA HiFi HotStart Real-Time PCR Master Mix is fully thawed and has been vortexed before use.

KAPA HiFi HotStart Real-Time PCR Master Mix (2X) may be stored at 4 °C for regular, short-term use (up to 1 month). Long-term storage at room temperature or 4 °C is not recommended. Please note that reagents stored above -20 °C are more prone to degradation when contaminated by the user; storage at such temperatures is therefore at the user's own risk.

### 6.2 Quality control

All kit components are subjected to stringent quality control tests, are free of contaminating exo- and endonuclease activities and meet strict requirements with respect to DNA contamination. Detailed quality control information for individual kit components is available upon request, please contact [support@kapabiosystems.com](mailto:support@kapabiosystems.com).

### 6.3 Product use limitations and licenses

KAPA Library Preparation Kits are developed, designed and sold exclusively for research purposes and *in vitro* use. Neither the product, nor any individual component, has been tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to the MSDS, which is available on request.

Certain applications of this product are covered by patents issued to parties other than Kapa Biosystems and applicable in certain countries. Purchase of this product does not include a license to perform any such applications. Users of this product may therefore be required to obtain a patent license depending upon the particular application and country in which the product is used.

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,994,056 and 6,171,785. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as apparatus or system claims in US Patent No. 6,814,934) and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

This product is provided under an agreement between Molecular Probes, Inc. and Kapa Biosystems Inc., and the manufacture, use, sale or import of this product is subject to one or more of U.S. Patent Nos. 5,436,134; 5,658,751 and corresponding international equivalents, owned by Molecular Probes, Inc. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer, where such research does not include testing, analysis or screening services for any third party in return for compensation on a per test basis. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851. SYBR® is a registered trademark of Molecular Probes, Inc.

The purchase of this product includes a limited, non-transferable license under specific claims of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, owned by the University of Utah Research Foundation or Evotec Biosystems GmbH and licensed to Idaho Technology, Inc. and Roche Diagnostics GmbH, to use only the enclosed amount of product according to the specified protocols. No right is conveyed, expressly, by implication, or by estoppel, to use any instrument or system under any claim of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, other than for the amount of product contained herein.

Licensed under U.S. Patent nos. 5,338,671 and 5,587,287 and corresponding patents in other countries. TruSeq is a trademark of Illumina, Inc. SYBR® is a registered trademark of Molecular Probes, Inc.

For technical support please contact [support@kapabiosystems.com](mailto:support@kapabiosystems.com)

#### Boston, Massachusetts, United States

600 West Cummings Park, Suite 2250  
Woburn, MA 01801 U.S.A.  
Tel: +1 781 497 2933 Fax: +1 781 497 2934

#### Cape Town, South Africa

2nd Floor, Old Warehouse Building, Black River Park,  
Fir Road, Observatory 7925, Cape Town, South Africa  
Tel: +27 21 448 8200 Fax: +27 21 448 6503

Email: [info@kapabiosystems.com](mailto:info@kapabiosystems.com)