



KAPA™ PROBE FAST qPCR Kit

Master Mix (2X) ABI Prism™

1. Product Description

KAPA PROBE FAST ABI Prism™ qPCR Kits are designed for high throughput, fast-cycling, real-time PCR using sequence-specific fluorogenic probes. These kits are compatible with all fluorogenic probe-based technologies, including hydrolysis probes (e.g. TaqMan®) and displacement probes (e.g. molecular beacons).

KAPA PROBE FAST ABI Prism™ qPCR Master Mix (2X) is a ready-to-use cocktail containing all components except primers, probe(s) and template for fast cycling probe-based real-time PCR on ABI real-time instruments that support normalization with ROX reference dye at a concentration of 500 nM. The 2X Master Mix contains KAPATaq HotStart DNA polymerase, KAPA PROBE FAST qPCR Buffer, dNTPs, MgCl₂, ROX reference dye and stabilizers.

KAPATaq HotStart DNA Polymerase is an antibody-mediated hot start formulation of KAPATaq DNA polymerase. In the HotStart formulation, the enzyme is combined with a proprietary antibody that inactivates the enzyme until the first denaturation step. This eliminates spurious amplification products resulting from non-specific priming events during reaction setup and initiation, and increases overall reaction efficiency.

2. Product Applications

KAPA PROBE FAST ABI Prism™ qPCR Kits are ideally suited for:

- Gene expression analysis
- SNP genotyping
- Microarray validation
- Gene knockdown validation

3. Product Specifications

3.1 Shipping and Storage

KAPA PROBE FAST ABI Prism™ qPCR Kits are shipped on ice packs. Upon arrival, store kit components protected from light at -20 °C in a constant-temperature freezer. When stored under these conditions and handled correctly, the 2X Master Mix will retain full activity for 18 months from the date of receipt.

3.2 Handling

ROX reference dye is sensitive to exposure to light. Avoid repeated freezing and thawing. Always ensure that the product has been fully thawed and mixed before use.

3.3 Quality Control

KAPA PROBE FAST ABI Prism™ qPCR Master Mix (2X) is free of contaminating DNase and RNase. It is functionally tested to demonstrate resolution of 5 orders of linear dynamic range using human genomic DNA as template and an ApoB100 primer/probe assay.

3.4 Product Use Limitations

KAPA PROBE FAST ABI Prism™ qPCR Master Mix (2X) is sold exclusively for research purposes and *in vitro* use. Neither the product, nor any individual component, was tested for use in diagnostic applications or for drug development, nor is it suitable for administration to humans or animals. Please refer to the MSDS, which is available upon request.

KK4705

 100 x 20 µl
reactions

KAPA PROBE FAST qPCR Master Mix (2X)
ABI Prism™ 1 x 1 ml

 Contains:
- qPCR Master Mix (2X) with ROX reference dye

KK4706

 500 x 20 µl
reactions

KAPA PROBE FAST qPCR Master Mix (2X)
ABI Prism™ 1 x 5 ml

 Contains:
- qPCR Master Mix (2X) with ROX reference dye

KK4707

 1000 x 20 µl
reactions

KAPA PROBE FAST qPCR Master Mix (2X)
ABI Prism™ 2 x 5 ml

 Contains:
- qPCR Master Mix (2X) with ROX reference dye

 The final MgCl₂ concentration per PCR reaction is 5 mM

Quick Notes

- This kit is designed for high-throughput, fast cycling, real-time PCR using sequence-specific fluorogenic probes.
- The kit is suitable for all fluorogenic probe-based technologies, including hydrolysis probes (e.g. TaqMan®) and displacement probes (e.g. molecular beacons).
- Initial denaturation for 20 sec at 95 °C is sufficient for enzyme reactivation, however optimal denaturation of complex targets may require up to 3 min denaturation.
- For two-step cycling, use 30 sec combined annealing/extension/data acquisition.
- For three-step cycling, use 20 sec for primer annealing and 1 sec for extension/data acquisition.
- Do not exceed 25 µl reaction volumes.

Instrument Table

Instrument	ROX reference dye
ABI 5700, 7000, 7300, 7700, 7900HT, StepOne™, StepOnePlus™ and ViiA™ 7	500 nM final



KAPA PROBE FAST qPCR Kit

Master Mix (2X) ABI Prism™

4. KAPA PROBE FAST qPCR Protocol

Any existing qPCR assay performed efficiently using standard cycling conditions may be converted to a Fast qPCR assay with KAPA PROBE FAST qPCR Kits. Typically, minimal re-optimization of reaction parameters is required.

This protocol is intended for use with the ABI PRISM™7000, 7700, 7900HT, the ABI 7300 Real-Time PCR System, the GeneAmp®5700, the StepOne™ and the StepOnePlus™ instruments. This kit is not compatible with instruments that use ROX at a final concentration lower than 500 nM.

4.1 Step 1: qPCR Reaction Setup

- Before preparing qPCR reactions, thoroughly mix the KAPA PROBE FAST ABI Prism™ qPCR Master Mix (2X), template DNA, primers and probes.
- Calculate the required volumes of each component based on the following table:

	Final concentration	20 µl rxn
PCR grade water up to 20 µl		As required
KAPA PROBE FAST qPCR Master Mix (2X) - ABI Prism™	1X	10 µl
Forward Primer (10 µM)	100 - 400 nM	Variable
Reverse Primer (10 µM)	100 - 400 nM	Variable
Probe	100 - 500 nM	Variable
Template DNA or cDNA	<250 ng	Variable

4.2 Step 2: Plate Setup

- Preparation of a reaction cocktail is vital in qPCR to reduce the effect of pipetting errors between samples. Assemble all components above except template DNA or cDNA.
- Gently mix all components in the cocktail before transferring the appropriate volume of reaction mixture to each well of a PCR tube/plate.
- Add template DNA or cDNA to each reaction.
- Reaction volumes may be scaled down from 20 µl to 10 µl if low volume tubes/plates are used.
- Cap or seal the reaction tube/plate and centrifuge briefly.

4.3 Step 3: Run the qPCR reaction

- If applicable, select fast mode on the instrument.
- Program the following cycling protocol:

Step	Temperature	Duration	Cycles
Enzyme activation	95 °C	20 sec - 3 min ¹	Hold
Denature	95 °C	1 - 3 sec	40
Anneal/Extend/Acquire	60 °C	≥20 sec ²	

¹20 sec at 95 °C is sufficient time for enzyme activation, however optimal denaturation of complex targets may require up to 3 min denaturation.

²Select minimum programmable time (not less than 20 sec) according to instrument guide.

4.4 Step 4: Analyze the results

- Data analysis varies depending on the instrument used. Please refer to your instrument user guide for information.



5. Important Parameters

5.1 Assay Design

We recommend using previously validated assays or using dedicated qPCR design software such as Beacon Designer 7 when designing Probe-based assays (www.PremierBiosoft.com).

Lyophilized primers and probes should be resuspended in 10 mM Tris-HCl pH 8.0, 1 mM EDTA. DNA kept frozen in a nuclease-free environment should be stable for years. We find it convenient to initially prepare a 100 μ M freezer stock (which should be thawed relatively infrequently).

Optimal primer concentration should be determined empirically. To maximize the sensitivity of the assay, use the lowest concentration of primers that can be used without compromising the efficiency of the qPCR reaction. The optimal primer concentration range is 100 – 400 nM.

Optimal probe concentration should be determined empirically. The optimal probe concentration range has generally been found to be 100 – 500 nM.

5.2 MgCl₂

The concentration of MgCl₂ affects the binding dynamics of primers and probes to template DNA. The higher the final MgCl₂ concentration in the PCR reaction, the greater the binding affinity of the primers and probe for target DNA. KAPA PROBE FAST ABI Prism™ qPCR Master Mix (2X) provides MgCl₂ at a final concentration of 5 mM, which is suitable for most targets.



6. Troubleshooting

Symptom	Possible Cause	Solution
Late C _t or no amplification during cycling	<p>Incorrect cycling protocol</p> <p>Incorrect reaction setup</p> <p>Incorrect detection filter/channel</p> <p>Degraded template DNA</p> <p>Sub-optimal primer/probe design</p> <p>Degraded probe or primers</p> <p>Non-specific products may be amplified</p> <p>Incorrect annealing temperature</p>	<p>Verify that the correct default cycling conditions were used.</p> <p>Verify that all the components have been added at the correct concentrations.</p> <p>Check that the correct filters have been selected for data acquisition.</p> <p>Prepare fresh template then repeat experiment.</p> <p>We recommend using pre-validated assays or designing them using dedicated software.</p> <p>Check integrity of primers/probe on a denaturing polyacrylamide gel.</p> <p>Repeat setup on ice and run qPCR reaction immediately after setup.</p> <p>Increase annealing temperature in 2 °C increments.</p> <p>Decrease annealing temperature in 2 °C increments.</p>
The NTC gives a positive result	<p>Contamination of reagents</p> <p>Contamination during setup</p> <p>Degradation of primers and probe</p>	<p>Discard all reagents and repeat experiment with new components.</p> <p>Review setup procedure and ensure that aerosol-barrier pipette tips are used.</p> <p>Use new stocks of primers and probe or redesign the assay.</p>
Extremely high ΔRn or Rn values	<p>ROX was not selected as the passive reference dye at setup</p>	<p>Select ROX as the passive reference dye when setting up the plate.</p>
High variability across replicates	<p>Insufficient mixing of reaction master mix</p> <p>Evaporation</p>	<p>Mix the reaction by inverting the tube a few times, followed by brief centrifugation prior to aliquotting to the reaction plate.</p> <p>Ensure that the optical lids or sealing film is completely sealed before loading qPCR instrument. This is particularly important on the edges of qPCR plates.</p>

7. Note to Purchaser: Limited License

Certain applications of this product are covered by patents issued to parties other than Kapa Biosystems and applicable in certain countries. Purchase of this product does not include a license to perform any such applications. Users of this product may therefore be required to obtain a patent license depending upon the particular application and country in which the product is used.

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Authorized 5' Nuclease Core Kit and Licensed Probe, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

This product is an Authorized 5' Nuclease Core Kit. Use of this product is covered by one or more of the following claims outside the US corresponding to US Patent No. 5,210,015, 5,487,972, and 5,928,907 (claim numbers 12-24, 27-28). The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. Separate purchase of a Licensed Probe would convey rights under US Patents and corresponding patent claims outside the US: 5,538,848, 5,723,591, 5,876,930, 6,030,787, 6,258,569, 5,804,375 (claims 1-12 only), and claims outside the United States corresponding to US Patent No. 6,214,979. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Licensed under U.S. Patent nos. 5,338,671 and 5,587,287 and corresponding patents in other countries.

Prism™, StepOne™, StepOnePlus™, ViiA™ 7, GeneAmp®, ROX and TaqMan® are trademarks or registered trademarks of their respective companies.

For technical support please contact:
support@kapabiosystems.com

Boston, Massachusetts, United States

600 West Cummings Park, Suite 5350
Woburn, MA, 01801 U.S.A.
Tel: +1 781 497 2933 Fax: +1 781 497 2934

Cape Town, South Africa

2nd Floor, Old Warehouse Building, Black River Park,
Fir Road, Observatory, 7925 Cape Town, South Africa
Tel: +27 21 448 8200 Fax: +27 21 448 6503

Email: info@kapabiosystems.com